Quality and Experience

Catheter Kit

(de) Katheter-Kit (fr) Kit de cathéter

(it) Kit catetere (es) Kit de catéteres

(en) Instructions for Use	2
(de) Gebrauchsanweisung	15
(fr) Mode d'emploi	31
(it) Istruzioni per l'uso	46
(es) Instrucciones de uso	60

www.pfmmedicalusa.com


PFM Medical, Inc.
1916 Palomar Oaks Way, Suite 150
Carlsbad, CA 92008
Tel.: (760) 758-8749
Fax: (760) 758-1167
www.pfmmedicalusa.com
customerservice@pfmmedicalusa1.com

pfm medical ag
Wankelstr. 60 50996 Köln,
Germany

Explanation of the symbols on label and packaging

	Do not re-use
	Do not use if package is damaged and consult instructions for use
	Use-by date
	Consult instructions for use or consult electronic instructions for use
	Keep away from sunlight
	Keep dry
	Non-pyrogenic
	Batch Code
	Do not re-sterilize
	Catalogue number
	Manufacturer
	Prescription only
	Sterilized using ethylene oxide
	Authorized representative in the European Community/ European Union
	Importer
	Distributor
	Medical Device
	Single sterile barrier system with protective packaging inside
	CE marking
	No Latex

	Patient Name or patient ID
	Date of implantation
	Name and Address of the implanting healthcare institution/provider
	Information website for patients
	Caution
	Unique Device Identifier

(en) Instructions for Use

***prime*PICC® Peripherally inserted central venous catheter**

Product Description

The *prime*PICC catheters are a family of peripherally inserted central catheters made from specially formulated biocompatible medical grade material. Each *prime*PICC catheter has a kink-resistant, reverse taper design. Catheters are packaged in a tray with accessories necessary for reliable long-term or short-term vascular access.

Indications for Use

3F Single Lumen, 4F Single Lumen, 5F Single Lumen, 5F Dual Lumen *prime*PICC Catheters:

The *prime*PICC catheter is indicated for short or long-term access to the central venous system for intravenous therapy, blood sampling and power injection of contrast media. The maximum recommended infusion rate is 5 ml/sec. The maximum pressure of power injectors used may not exceed 300 psi.

4F Dual Lumen *prime*PICC Catheter:

The 4F Dual Lumen *prime*PICC catheter is indicated for short or long-term access to the central system for intravenous therapy, blood sampling and power injection of contrast media. The maximum recommended infusion rate is 4 ml/sec. The maximum pressure of power injectors used with the 4F Dual Lumen *prime*PICC catheter may not exceed 300 psi.

5F and 6F Triple Lumen *prime*PICC Catheters:

The *prime*PICC catheter is indicated for short or long-term access to the central venous system for intravenous therapy, blood sampling and power injection of contrast media. The maximum recommended infusion rate is 5 ml/sec. The maximum recommended infusion rate is 5 ml/sec. The maximum pressure of power injectors used may not exceed 300 psi.

The primary line of the 5F and 6F *prime*PICC Catheter is also indicated for central venous pressure monitoring. For central venous pressure monitoring, it is recommended that a catheter of lumen 20 gauge or larger be used.

Contraindications

The *prime*PICC is contraindicated whenever:

- The patient's body size or vein size is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device. The implantable portion of the catheter is composed of 80% polyurethane and 20% barium sulfate.
- There has been past irradiation of prospective insertion site.
- There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.
- There are local tissue factors that may prevent proper device stabilization and/or access. Or skin conditions that prevent insertion through an unaffected area.
- There is systemic infection or local infection at the insertion site.

General Information and Warnings

The product described here may be used only by specialist physicians who are familiar with the techniques used in diagnostic and interventional cardiology. No additional training by PFM Medical is required for use of the *prime*PICC.

These instructions for use and the information on the packaging should be read carefully before each use.

Included with the product is an Implant Card to be filled out and given to the patient containing the product information. Please fill in the following information prior to giving to the patient: name of the patient or patient ID, date of implantation, and name and address of healthcare intuition/provider.

The target treatment population for the *prime*PICC are patients who are clinically deemed appropriate to accommodate the device size including neonates, children, adolescents, and adults.

The Summary of Safety and Clinical Performance is available on Eudamed.

The *prime*PICC products have a maximum shelf life of 3 years.

General Warnings

- This catheter is for Single Use Only. Do not reuse. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.
- Do not re-sterilize the catheter or accessories by any method.
- The manufacturer shall not be liable for any damages caused by reuse or re-sterilization of this catheter or accessories.
- Contents sterile and non-pyrogenic in unopened, undamaged package.
- Sterilized by ethylene oxide.
- Do not use catheter or accessories if package is opened or damaged.
- Do not use catheter or accessories if any sign of product damage is visible.
- Catheter must not be implanted in infectious area and in patients with systemic infection.
- The *prime*PICC catheter features a reverse-taper catheter design. Placement of a larger catheter at or below the antecubital fossa may result in an increased incidence of phlebitis. Placement of the *prime*PICC catheter above the antecubital fossa is recommended.
- When connecting accessories to the *prime*PICC catheters, only use devices with an ISO compliant luer.

Placement Warnings

- In the rare event that a hub or connector separates from any component during insertion or use, take all necessary steps and precautions to prevent blood loss or air embolism and remove the catheter.
- Do not advance the guidewire or catheter if unusual resistance is encountered.
- Do not insert or withdraw the guidewire quickly or forcibly from any component. The wire may break or unravel. If the guidewire becomes damaged, the introducer needle or sheath/dilator and guidewire must be removed together.

Power Injection Warnings

- Failure to warm contrast to body temperature prior to power injection may result in catheter failure.
- Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure.

- Use of lumens not marked “Power Injectable” for power injection of contrast media may cause failure of the catheter.
- Power injector machine pressure limiting feature may not prevent over pressurization of an occluded catheter.
- Exceeding the maximum indicated flow rate may result in catheter failure and/or catheter tip displacement.
- *prime*PICC catheter indication of power injection of contrast media implies the catheter's ability to withstand the procedure but does not imply appropriateness of the procedure for a particular patient. A suitably trained clinician is responsible for evaluating the health status of a patient as it pertains to a power injection procedure.

Precautions

- Small syringes will generate excessive pressure and may damage the catheter. Ten (10) cc or larger syringes are recommended.
- Do not use sharp instruments near the extension lines or catheter lumen.
- Do not use scissors to remove dressing.
- Catheter will be damaged if clamps other than what is provided with this kit are used.
- Clamping of the tubing repeatedly in the same location will weaken tubing. Avoid clamping near the luer(s) and hub of the catheter.
- Examine catheter lumen and extension(s) before and after each infusion for damage.
- To prevent accidents, assure the security of all caps and connections prior to and between uses.
- Use only Luer Lock (threaded) Connectors with this catheter.
- Repeated over tightening of luer lock connections, syringes, and caps will reduce connector life and could lead to potential connector failure.
- Confirm catheter tip position by x-ray prior to use. Monitor tip placement routinely per institution policy.

Possible Complications

Insertion of these devices requires training consistent with the Standards of Practice. Ensure that you are familiar with the below complications and their emergency treatment should any of them occur.

- | | |
|-------------------------------------|--|
| • Air Embolism | • Hematoma |
| • Accidental Arterial Access | • Hemolyzed blood |
| • Bleeding | • Intolerance Reaction to Implanted Device |
| • Brachial Plexus Injury | • Laceration/Perforation of Vessels or Viscus |
| • Cardiac Arrhythmia | • Myocardial Erosion |
| • Cardiac Tamponade | • Phlebitis |
| • Catheter Breakage, leak, burst | • Spontaneous Catheter Tip Malposition or Retraction |
| • Catheter Erosion through the Skin | • Thromboembolism |
| • Catheter Related Sepsis | • Venous Thrombosis |
| • Endocarditis | • Ventricular Thrombosis |
| • Exit Site Infection | • Risks Normally Associated with local or General Anesthesia, Surgery, and Post-Operative Recovery |
| • Exit Site Necrosis | |
| • Extravasation | |
| • Fibrin Sheath Formation | |

Directions for Use

1. Identify the Vein and Insertion Site

- a. Open insertion kit and access kit/supplies using sterile technique.
- b. Apply tourniquet above anticipated insertion site.
- c. Select vein based on patient assessment. Recommended veins are the basilic, cephalic and median cubital veins.
- d. Release tourniquet.

2. Patient Position / Catheter Measurement

- a. Position arm at a 90° angle.
- b. Measure from the planned insertion site to the right clavicular head, then down to the third intercostal space.

Note: For central venous placement, the recommended target tip location is the lower 1/3 of the Superior Vena Cava or Cavoatrial Junction. The external measurement does not exactly duplicate the internal venous anatomy.

3. Skin Preparation

- a. Don sterile gloves.
- b. Apply under drape.
- c. Prepare the site with included skin antiseptic, using sterile technique.
- d. If alcohol is used, allow to air dry completely before proceeding with procedure.
- e. Remove/discard gloves.

4. Sterile Field Preparation

- a. Apply tourniquet above the intended insertion site.
- b. Don sterile gloves.
- c. Apply drapes and complete sterile field preparation using sterile technique.

5. Prepare Catheter and Hydrate Stylet

- a. Attach prefilled syringe to the luer attachment on the Touhy Borst attachment (if included) or to lumen of catheter and flush catheter with sterile saline.
- b. If using double lumen catheter, attach needleless access port to remaining extension. Attach saline filled syringe to the needleless access port and completely flush catheter lumen. Remove syringe from needleless access port prior to clamping extension.
- c. The syringe may be left attached during the procedure.

Caution: Never close clamp on catheter stylet; stylet and catheter damage may result.

Caution: The needleless access port should not be used with needles, blunt cannula, or other non-luer connectors, or luer connectors with visible defects. If needle access is attempted, the needleless access port must be replaced immediately. Do not exceed 100 actuations.

6. Modification of Catheter Length

- a. Measure the distance from the insertion site (zero mark) to the desired tip location, replicating measurement as described in step 2.

Note: If placing *prime*PICC using fluoroscopic guidance, this step will occur after guidewire is inserted through introducer needle and advanced to intended tip location (step 7 h).

- b. Disconnect Touhy Borst Adapter from luer hub (if included).
Note: If catheter kit does not contain the Touhy Borst Adapter and stylet, skip to step 7.
- c. Withdraw entire Touhy Borst Adapter assembly as one unit.
- d. Retract the stylet to well behind the point that the catheter is to be cut.
- e. Using sterile scalpel or scissors, carefully cut the catheter according to institutional policy if necessary.
- f. Inspect cut surface to assure there is no loose material.
- g. Re-advance the Touhy Borst Adapter/Stylet to the luer hub. Assure stylet is intact.
Note: If catheter kit does not contain the Touhy Borst Adapter and stylet, skip to step 7.
- h. Gently retract the stylet through the Touhy Borst Adapter until the stylet is contained inside the catheter.
- i. Assure proper alignment of the stylet to the distal end of the trimmed catheter.
Caution: Never attempt to cut stylet.
Caution: Always withdraw stylet back beyond the tip of the catheter prior to insertion.

7. Perform Venipuncture

- a. Anesthetize the proposed insertion site with local anesthesia per institutional policy.
- b. Insert the introducer needle with attached syringe into the target vein and observe for blood flashback.
- c. Release tourniquet.
- d. Remove syringe from introducer needle and place thumb over needle hub to prevent air embolism.
- e. Insert the flexible end of marked guidewire with forward motion into and past the needle hub into the vein.
- f. Advance guidewire to desired depth.
Note: Do not advance the wire past the axilla without fluoroscopic guidance or other tip locating methods.
Caution: The length of the wire inserted is determined by the size of the patient and should not be advanced past the shoulder except when fluoroscopic guidance is available. Monitor patient for arrhythmia throughout this procedure. The patient should be placed on a cardiac monitor during this procedure. Cardiac arrhythmias may result if guidewire is allowed to pass into the right atrium. The guidewire should be held securely during this procedure.
- g. Remove needle, leaving guidewire in the target vein. Hold guidewire in position as needle is removed.
- h. Advance the micro introducer assembly over the guidewire by using a twisting motion. Continue advancement of micro introducer assembly into the vessel. If necessary, a small incision may be made adjacent to the guidewire to facilitate insertion of the sheath and dilator. Verify institutional guidelines regarding using scalpel for these purposes prior to attempting additional incision.
Caution: DO NOT bend the sheath/dilator during insertion as bending will cause the sheath to prematurely tear. Hold sheath/dilator close to the tip (approximately 3cm from tip) when initially inserting through the skin surface. To progress the sheath/dilator towards the vein, regrasp the sheath/dilator a few centimeters (approximately 5cm) above the original grasp location and push down on the sheath/dilator. If necessary, nick the skin with the safety scalpel for easier advancement of the dilator.

Note: If placing *prime*PICC using fluoroscopic guidance, advance the guidewire within the micro introducer until it reaches the lower 1/3 of the Superior Vena Cava or Cavoatrial Junction and confirm final tip position. Once final tip position is confirmed, clamp the guidewire with forceps and measure the depth of the guidewire by reading the markings on the guidewire. Modify the catheter length based on the guidewire measurement.

- i. Withdraw the dilator and guidewire, leaving the sheath in place.

8. Insert and Advance Catheter

- a. Place a finger over the orifice of the sheath to help prevent air embolism and excessive blood loss. Consider having patient perform Valsalva maneuver until catheter is inserted into sheath.
- b. Advance catheter slowly into and through the sheath.
- c. Continue to advance catheter. For central placement, when the catheter has advanced to the shoulder, have the patient turn head and place chin on shoulder toward the insertion site to prevent possible insertion into jugular vein.

Note: If placing *prime*PICC using fluoroscopic guidance and over-the-wire approach, insert proximal end of guidewire into distal end of catheter, and feed catheter over the guidewire into the vessel. Make sure the guidewire extends out of the proximal hub of the catheter before advancing the guidewire into the vessel. Advance catheter over the guidewire until distal tip is correctly positioned. Maintain control of guidewire external to proximal hub of catheter at all times.

- d. Due to the reverse taper on the proximal end of the catheter, you will feel resistance in the sheath approximately 3cm distal to the catheter hub. The introducer may be partially split, but not removed, to facilitate insertion of the catheter past this point if necessary.
- e. Complete catheter advancement into the desired position.

Note: If placing *prime*PICC using fluoroscopic guidance with over-the-wire approach, remove the wire from the catheter by applying gentle pressure with one hand above the insertion site while grasping the wire with the other hand and pulling gently back with a constant motion.

9. Removing Introducer Sheath

- a. Stabilize the catheter position by applying pressure to the vein distal to the introducer sheath.
- b. Withdraw the introducer sheath from the vein and away from the insertion site.
- c. Split the introducer sheath and peel it away from the catheter.

Caution: Do not pull apart the portion of the sheath that remains in the vessel. To avoid vessel damage and skin laceration, pull back the sheath as far as possible and tear the sheath only a few centimeters at a time. Never leave sheath in place as an indwelling catheter. Damage to the vein will occur.

10. Remove the Touhy Borst Adapter Assembly (when included with catheter)

- a. Disconnect the Touhy Borst Adapter Assembly from the catheter luer connector.
- b. Stabilize the catheter position by applying light pressure to the vein distal to the insertion site.
- c. Gently remove the Adapter Assembly and stylet as one unit. Do not remove stylet through Touhy Borst Adapter Assembly.
- d. Clamp catheter.

Caution: To avoid difficulty and/or bunching of the catheter lumens while removing the stylet, flush catheter prior to removal. The catheter may need to be repositioned to allow for removal of the stylet.

Caution: Do not attempt to reinsert stylet once it has been withdrawn.

Caution: Never leave stylet in place after catheter insertion; injury may occur. Remove stylet after insertion.

11. Aspirate and Flush Catheter

- a. Attach saline-filled syringe to hub.
- b. Aspirate for adequate blood return and flush catheter with 10ml sterile normal saline to ensure patency. Clamp catheter prior to disconnecting saline to prevent air embolism.
- c. Attach cap to hub

NOTE: Included MaxPlus Clear needle-free connector may also be attached. Refer to step 12 or MaxPlus Clear IFU for use.

Caution: Small syringes will generate excessive pressure and may damage the catheter. Ten (10) cc or larger syringes are recommended.

12. Use of MaxPlus Clear Needle-Free Connector

Indications for Use:

The MaxPlus is a sterile single patient use positive displacement connector for needleless access to the IV line and/or IV catheter during IV therapy. The MaxPlus connector can be used for direct injection, intermittent infusion, continuous infusion or aspiration.

Description:

The MaxPlus needleless connector is a closed luer activated device. The accessing ISO male luer from standard administration sets, extension sets, and syringes activate the flow of fluid through the device. The MaxPlus features Tru-Swab® technology which provides a flat, smooth surface for optimum disinfection during pre-access swabbing. The positive displacement feature of the MaxPlus product produces a positive bolus of fluid to clear the catheter upon disconnection of the male luer.

- a. Without touching the blue top remove the MaxPlus from the package. Do not contaminate.
- b. Attach primed IV set or syringe with flush/priming solution to the MaxPlus. Invert the MaxPlus and prime in accordance with facility protocol.
- c. Remove the protective cover and attach MaxPlus to the hub of the extension set or vascular access device by inserting luer into the hub and rotating until secure. Do not over tighten.
- d. Prior to every access, always swab the top of the MaxPlus with an appropriate antiseptic and allow to dry.
- e. Attach luer from primed IV set or syringe to the MaxPlus. If luer is a two piece spin collar, first pull back the collar and insert luer in a straight in motion and rotate ¼ turn clockwise; then push the spin collar forward and tighten. If syringe has a luer slip, insert and rotate ¼ turn clockwise to secure connection. Do not leave luer slip unattended.
- f. Flush the MaxPlus after each use with normal saline or in accordance with facility protocol.
- g. To disconnect male luer from MaxPlus, hold MaxPlus and rotate mating luer until disconnected. Upon removal of the luer from the MaxPlus a positive displacement of fluid will occur to reduce retrograde flow into the catheter lumen. If blood or fluid residuals remain on the connector surface after access, swab surface.
- h. Apply clamp, after disconnecting luer.

Warnings and Cautions:

- The device should not be used with needles, blunt cannula systems, non ISO luer connections, or luer connections with visible defects. Doing so may result in leakage and/or failure of the device.

- Luer slip connections should not be left unattended due to potential for disconnection.
- For proper use, clinicians must be familiar with and trained in the use of the MaxPlus device. Its use should be preceded by an established facility protocol. The device should be disinfected with an appropriate antiseptic agent, such as 70% IPA, prior to each access.
- The device is intended for use with ISO luer lock and luer slip connectors provided on standard IV administration sets, extension sets, and syringes.
- Trace lines before connection. Verify the line being connected to the MaxPlus is the appropriate Intravenous Therapy line.
- Failure to properly prime the device can result in reflux.
- The MaxPlus should be changed according to facility protocol or in accordance with currently recognized guidelines for IV therapy.
- Sterile disposable single patient use device may be accessed multiple times per hospital protocol. Reuse, reprocessing or re-sterilizing may lead to patient infection or other illness/injury.
- Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.

13. Verify Placement

- PICCs should be positioned with the catheter tip in the lower 1/3 of the Superior Vena Cava or Cavoatrial Junction. Verify correct catheter tip position using radiography or other appropriate technology.

Caution: Do not clamp the lumen portion of the catheter. Clamp only the extension(s). Do not use the serrated forceps; use only the in-line clamp(s) provided. Caution: Failure to verify catheter placement may result in serious trauma or fatal complication.

Note: If there is no blood return, correct the catheter position before use.

14. Catheter Securement and Wound Dressing

The StatLock® stabilization device is included in *primePICC®* catheter insertion kits. Please refer to the below instructions or Instructions for Use on the proper use and removal of the device. The StatLock® catheter stabilization device should be monitored daily and replaced every seven days.

StatLock® Safety and Efficacy Considerations

The StatLock® stabilization device is for single use only. Rx only. Sterile only when included as a component in a sterile kit. This package does not maintain sterility; use once catheter kit package is opened. Discard if package is opened or damaged.

Do not alter the StatLock® stabilization device or components. Procedure must be performed by trained personnel with knowledge of anatomical landmarks, safe technique, and potential complications.

StatLock® Contraindications

- Known tape or adhesive allergies

StatLock® Warning and Precautions

- Do not use the StatLock® stabilization device where loss of adherence could occur, such as with a confused patient or nonadherent skin.
- Observe universal blood and body fluid precautions and infection control procedures, during application and removal of the StatLock® stabilization device.
- Avoid the StatLock® stabilization device contact with alcohol or acetone, both can weaken bonding of components and the StatLock® stabilization device pad adherence.

- Minimize catheter/tube manipulation during application and removal of the StatLock® stabilization device.
 - Always apply adhesive strip to central venous catheters at or near insertion site during dressing change involving the StatLock® stabilization device.
 - Remove oil and moisturizer from targeted skin area.
 - The StatLock® stabilization device should be monitored and replaced in accordance with facility protocol.
 - **Warning:** Intended for single use. Do not reuse. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness, or death of the patient.
 - **Warning:** Do not re-sterilize. The sterility of the single use device is not guaranteed following re-sterilization because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Re-sterilization may compromise the structural integrity, essential material and/or design characteristics and may lead to an unpredictable loss of functionality and/or device failure.
 - Federal (USA) law restricts this device to sale by or on the order of a physician.
- a. Prepare target area with skin prep and allow to completely dry.
 - b. Place suture holes in catheter wing over StatLock® stabilization device posts (prior to pad adherence). Close lids, peel off paper backing, and place on target skin area.
 - c. To remove, lift plastic lids one side at a time.
 - d. Remove the StatLock® stabilization device anchor pad adhesive with alcohol.

Note: The StatLock® Stabilization device pad can be adhered after confirmation of proper catheter/tube position.

Note: Apply adhesive strip to central venous catheters for added security during dressing changes involving the StatLock® stabilization device.

Warning: The insertion site and external portion of the catheter should always be covered with a protective dressing.

15. Power Injection Procedure

- a. Remove the injection/needleless cap from the *prime*PICC catheter.

Note: If using included MaxPlus Clear injection/needleless cap, refer to step 12 or device IFU prior to removal.
- b. Attach a 10ml or larger syringe filled with sterile normal saline. Flush to ensure patency.

Note: Resistance to flushing may indicate partial or complete catheter occlusion. Do not proceed with power injection study until occlusion has been cleared.

Note: Vigorously flush the *prime*PICC catheter using a 10cc or larger syringe and sterile normal saline prior to and immediately following the completion of power injection studies. This will ensure the patency of the catheter and prevent damage to the catheter.
- c. Disconnect Syringe and clamp catheter.

Warning: Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure.
- d. Attach the power injection device to the *prime*PICC catheter per manufacturer's recommendations.
- e. Contrast media should be warmed to body temperature prior to power injection.

Warning: Always use injectable luer locking tubing between power injector syringe and catheter. Do not attempt to connect power injector syringe directly to the catheter. Damage may result.

- f. Use only lumens marked “Power Injectable” for power injection of contrast media.

Warning: Use of lumens not marked “Power Injectable” for power injection of contrast media may cause failure of the catheter.

- g. Unclamp catheter prior to beginning injection of contrast media.
- h. Complete power injection study per institution protocol taking care not to exceed the flow rate printed on the catheter of 300 psi at 5 ml/sec (for all configurations except for 4 French Dual Lumen) or 300 psi at 4 ml/sec (4 French Dual Lumen configurations).

Warning: Exceeding the maximum indicated flow rate may result in catheter failure and/or catheter tip displacement.

- i. Disconnect the power injection device.
- j. Attach a 10ml or larger syringe filled with sterile normal saline. Flush to ensure patency.
- k. Disconnect syringe and clamp catheter.

16. Dressing Changes/Exit Site Cleaning

- a. Assess the dressing for the accumulation of blood, fluid or moisture beneath the dressing in the first 24 hours. Assess the external length of the catheter to determine if migration of the catheter has occurred. Confirm catheter placement, tip location, patency, and security of dressing periodically.

Note: Scissors should not be used to remove the dressing changes to minimize the risk of cutting the catheter. Do not suture through or around any part of the catheter shaft or extension legs. If using sutures to secure the catheter, use the suture wings.

Note: When using alcohol or alcohol containing antiseptics with the *prime*PICC catheter, care should be taken to avoid prolonged or excessive contact. Solutions should be allowed to completely dry before applying an occlusive dressing. Chlorhexidine gluconate and/or povidone iodine are the suggested antiseptics to use.

Warning: Alcohol should not be used to soak or declog the *prime*PICC catheter because alcohol is known to degrade polyurethane catheters over time with repeated and prolonged exposure.

Warning: Acetone and polyethylene glycol containing ointments should not be used with the *prime*PICC catheter, as these may cause failure of the device.

18. Flushing and Heparinization

- a. Follow institutional policy for flushing frequency and heparin utilization.
- b. After drug administration each lumen should be flushed again with normal saline and then locked per institutional policy to maintain patency.
- c. Blood Sampling/TPN Infusions – follow institutional policy for flushing and locking procedure.
- d. Injection cap(s) or needleless access port(s) should be changed per institutional policy.

19. Occluded/Partially Occluded Catheter

- a. If resistance is encountered to aspirating or flushing, the lumen may be partially or completely occluded. If the lumen will neither aspirate nor flush, and it has been determined that the catheter is occluded, follow institutional de clotting procedure.

Warning: Do not flush against resistance.

20. Central Venous Pressure Monitoring (5 French and 6 French Triple Lumen Configurations)

- a. Prior to central venous pressure monitoring:
- b. Ensure catheter tip is properly positioned.
- c. Vigorously flush catheter with sterile normal saline.
- d. Ensure pressure transducer is at the level of the right atrium.
- e. A continuous infusion of saline (3ml/hr) is recommended while measuring CVP to improve accuracy of results
- f. Use institutional protocol for central venous pressure monitoring procedures.

Note: CVP monitoring should be used in conjunction with other patient assessment metrics when evaluating cardiac function.

21. Catheter Removal

Warning: Only a clinician familiar with the appropriate techniques should attempt the following procedures.

Caution: Always review facility protocol, potential complications and their treatment, warnings, and precautions prior to catheter removal.

- a. Remove dressing and StatLock® catheter stabilization device or tape securement.
Note: Do not use scissors to remove dressing to minimize the cutting of the catheter.
- b. Grasp catheter near insertion site and using a slow steady motion, remove catheter from vein. Do not use excessive force.
- c. If resistance is felt – STOP. Apply a warm compress to the extremity for 20-30 minutes and refer to institutional policy.
- d. If no resistance is felt- Resume removal procedure.
- e. Apply pressure, if necessary, until bleeding stops and dress site following institutional policy.

Note: Inspect catheter and measure length, to ensure it matches documented length at insertion.

22. Disposal after Use

After use, medical products and accessories pose a potential biological hazard. For this reason, the products and their accessories should be handled and disposed of in accordance with recognized medical procedure, and in compliance with the relevant legal regulations and local ordinances.

Notice

Report any serious incident that has occurred in relation to the device to PFM Medical, Inc. and the competent authority of the Member State in which the user and/or patient is established.

Warranty

PFM Medical, Inc. warrants that this product was manufactured according to applicable standards and specifications. Patient condition, clinical treatment, and product maintenance may affect the performance of this product. Use of this product should be in accordance with the instructions provided and as directed by the prescribing physician.

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