Quality and Experience

Catheter Kit

(de) Katheter-Kit (fr) Kit de cathéter

(it) Kit catetere (es) Kit de catéteres

(en) Instructions for Use	2
(de) Gebrauchsanweisung	23
(fr) Mode d'emploi	47
(it) Istruzioni per l'uso	69
(es) Instrucciones de uso	91

k

www.pfmmedicalusa.com

PFM Medical, Inc.
1916 Palomar Oaks Way, Suite 150
Carlsbad, CA 92008
Tel.: (760) 758-8749
Fax: (760) 758-1167
www.pfmmedicalusa.com
customerservice@pfmmedicalusa1.com

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Explanation of the symbols on label and packaging

	Do not re-use
	Do not use if package is damaged and consult instructions for use
	Use-by date
	Consult instructions for use or consult electronic instructions for use
	Keep away from sunlight
	Keep dry
	Non-pyrogenic
	Batch Code
	Do not re-sterilize
	Catalogue number
	Manufacturer
	Date of Manufacture
	Sterilized using ethylene oxide
	Authorized representative in the European Community/ European Union
	Importer
	Medical Device
	Single sterile barrier system with protective packaging inside
	CE marking
	No Latex
	Patient Name or patient ID

	Date of implantation
	Name and Address of the implanting healthcare institution/provider
	Information website for patients
	Caution
	Unique Device Identifier
	MR Safe

(en) Instructions for Use

*prime*MIDLINE® catheter

Product Description

The *prime*MIDLINE Catheters are a family of peripherally placed catheters made from soft polyurethane materials. Each *prime*MIDLINE Catheter is packaged in a tray with the accessories necessary to establish short-term (less than 30 days) vascular access. The *prime*MIDLINE Catheters are suitable for use with power injectors.

Indications for Use

The *prime*MIDLINE Catheters are indicated for short term access to the peripheral venous system for clinician selected intravenous therapies, blood sampling and power injection of contrast media. These catheters are indicated for patients older than 28 days and that weigh more than 10 kg, with consideration given to adequacy of vascular anatomy and appropriateness of the procedure. The *prime*MIDLINE Catheters are suitable for use with power injectors.

Contraindications

This device is contraindicated whenever:

- The patient is under 28 days of age or under 10 kg in weight.
- The patient's body size or vein size is insufficient to accommodate the size of the implanted device.
- Bacteremia, septicemia, or other device-related infections are present or suspected.
- The patient is known or is suspected to be allergic to materials contained in the device.
- There has been past irradiation of prospective insertion site.
- There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement/body site or on the same extremity.
- There are local tissue factors that may prevent proper device stabilization and/or access, or skin conditions that prevent insertion through an unaffected area.

General Information and Warnings

The product described in this Instructions for Use may be used only by clinicians who are familiar with the techniques and clinical practices. The Instructions for Use and the information on the packaging should be read carefully before each use.

General Warnings

- Intended for Single Use. DO NOT REUSE. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure and/or lead to injury, illness or death of the patient.
- Sterilized by Ethylene Oxide. Do not re-sterilize the catheter or accessories by any method.
- When using alcohol or alcohol-containing antiseptics with polyurethane *prime*MIDLINE Catheters, care should be taken to minimize contact duration. Solutions should be allowed to completely dry before applying an occlusive dressing. Use all antiseptics per hospital protocol.
- Acetone-based solutions, tincture of iodine or polyethylene glycol-containing ointments should not be used with polyurethane catheters as they may cause failure of the device.
- Alcohol should not be used to lock, soak, or de-clot polyurethane *prime*MIDLINE Catheters because alcohol is known to degrade polyurethane catheters over time with repeated and prolonged exposure.

- Contents are sterile and non-pyrogenic in an unopened, undamaged package. Do not use catheter or accessories if package is opened or damaged or if any sign of product damage is visible.
- In the rare event that a hub or connector separates from any component during insertion or use, take all necessary steps and precautions to prevent blood loss or air embolism and remove the catheter.
- Do not use the catheter if there is any evidence of mechanical damage or leaking. Damage to the catheter may lead to rupture, fragmentation, possible embolism, and surgical removal.
- If signs of extravasation exist, discontinue injections. Begin appropriate medical intervention immediately.
- The fluid level in the catheter will drop if the catheter connector is held above the level of the patient's heart and opened to air.

Placement Warnings

- Do not advance the guidewire or catheter if unusual resistance is encountered. Do not advance the guidewire past the axilla.
- Do not insert or withdraw the guidewire quickly or forcibly from any component. If the guidewire must be withdrawn while the needle is inserted, remove both needle and wire as a unit to prevent the needle from damaging or shearing the guidewire.
- Do not use excessive force when introducing the microintroducer as this can lead to vessel perforation, bleeding, and an increase risk for thrombosis occurrence.
- Ensure the stylet tip does not extend beyond the trimmed end of the catheter. Extension of the stylet tip beyond the catheter end may result in vessel damage, stylet damage, difficult removal, stylet tip separation, potential embolism and/or risk of patient injury. Do not cut the stylet.
- Placement of the *prime*MIDLINE Catheter above the antecubital fossa is recommended.
- (Pediatric) Insertion techniques and placement locations are often modified according to the size and developmental age of the child. Only clinicians experienced in proper positioning and placement of midlines in pediatric patients should place the catheter in this patient population.
- (Pediatric) This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion or cardiac tamponade.

Warnings Specific to Power Injection

- Contrast media should be warmed to body temperature prior to power injection.
- Failure to warm contrast to body temperature prior to power injection may result in catheter failure, vessel damage, and foreign body embolism.
- Vigorously flush the *prime*MIDLINE Catheter using a 10 mL or larger syringe and sterile normal saline (0.9% sodium chloride solution) prior to and immediately following the completion of power injection studies. This will ensure the patency of the catheter and prevent damage to the catheter. Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure, vessel damage, and foreign body embolism.
- Resistance to flushing may indicate partial or complete catheter occlusion. Do not proceed with power injection study until occlusion has been cleared.
- The pressure limiting feature of the power injector machine may not prevent over pressurization of an occluded catheter.
- Exceeding the maximum flow rate or the maximum pressure of power injectors of 300 psi may result in catheter failure, vessel damage, and foreign body embolism.
- The *prime*MIDLINE Catheter indication for power injection of contrast media implies the catheter's ability to withstand the procedure but does not imply the appropriateness of the procedure for a particular patient. A suitably trained clinician is responsible for evaluating the health status of a patient as it pertains to a power injection procedure.

Precautions

- Follow all contraindications, warnings, cautions, precautions and instructions for all infusates as specified by their manufacturer.
- The selection and appropriateness of intravenous therapies is the responsibility of the health care practitioner. The toxicity of the infusates including cytotoxic vesicant drugs shall be considered.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with medical practice and applicable local, state and federal laws and regulations.
- Only qualified health care practitioners should insert, manipulate and remove these devices.
- Use aseptic techniques whenever the catheter lumen is opened or connected to other devices.

Precautions Specific to Device Placement Procedure:

- Consider an alternate placement site when there has been past irradiation of the prospective insertion site or previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.
- Measure the catheter-to-vessel ratio prior to insertion of an upper extremity Vascular Access Device (VAD); ensure a catheter-to-vessel ratio of less than 45% [INS, 2021].
- Small syringes will generate excessive pressure and may damage the catheter. Do not use a syringe smaller than 10 mL to flush and confirm patency. Patency should be assessed with a 10 mL or larger syringe with sterile saline. Follow facility guidelines for flushing volumes. Do not infuse against resistance.
- Do not cut guidewire to alter length.
- Do not insert the stiff end of guidewire into the vessel as this may result in vessel damage.
- Keep sufficient guidewire length exposed at the hub to allow for proper handling. A non-controlled guidewire can lead to wire embolism.
- Never leave the stylet or stiffening wire in place after catheter insertion.
- Never use force to remove the stylet. Resistance can damage the catheter. If resistance or bunching of the catheter is observed, stop stylet withdrawal and allow the catheter to return to normal shape. Withdraw both the catheter and stylet together approximately 2 cm and reattempt stylet removal. Repeat this procedure until the stylet is easily removed. Once the stylet is out, advance the catheter to the desired position (zero mark).
- Catheter tip should terminate at the level of the axilla in children and adults [INS, 2021].
- Do not clamp extension leg when stylet or stiffening wire is in the catheter to minimize the risk of catheter or component damage.
- Do not use sharp instruments near the extension lines or catheter lumen.
- Do not suture through or around any part of the catheter's tubing (shaft or extension legs). If using sutures to secure the catheter, use the suture wings and make sure the sutures do not occlude, puncture or cut the catheter.
- Catheter will be damaged if clamps other than what is provided with this kit are used.
- Clamping of the tubing repeatedly in the same location will weaken tubing. Avoid clamping near the luer(s) and hub of the catheter.
- Examine catheter lumen and extension(s) before and after each infusion for damage.
- To prevent accidents, assure the security of all caps and connections prior to and between uses.
- Use only luer lock (threaded) connectors with this catheter.
- Repeated overtightening of luer lock connections, syringes, and caps will reduce connector life and could lead to potential connector failure.
- The catheter must be secured in place to minimize risk of catheter breakage and embolization. Do not use scissors to remove dressing.

Special Patient Population

NOTE: (Pediatric) "Provide vascular access with attention to the child's anatomical, physiological, and developmental level." [INS, 2021].

NOTE: (Pediatric) Midline catheters are inserted into a peripheral vein of the upper arm via the basilic, cephalic, or brachial vein with the terminal tip located at the level of the axilla in children and adults. [INS, 2021].

NOTE: (Pediatric) Follow manufacturer's recommendations and institution's policies, procedures and/or practice guidelines when prepping the insertion site and surrounding skin; and in the use of any drugs or medications such as chlorhexidine prep solutions, lidocaine injections and heparin flush solutions. Insertion of midlines in pediatric patients may require the use of accessories or components not included in this kit configuration, based on the size and developmental age of the child and facility protocol.

NOTE: (Pediatric) "The type and duration of infusion therapy, patient preference, and the patient's physiologic condition (eg, age, diagnosis, comorbidities) and vascular condition (eg, history of vascular access attempts, vessel and skin health at site of insertion and proximal) are assessed when preparing for site selection and VAD insertion." [INS 2021]. In addition, facility policies, procedures, and/or practice guidelines can be used to access proper site selection.

NOTE: (Pediatric) Prep the insertion site and surrounding skin per facility policies, procedures and/or practice guidelines.

Undesirable Side-Effects

Insertion of these devices requires training consistent with the Standards of Practice. Ensure that you are familiar with possible undesirable side-effects and their emergency treatment should any of them occur.

- Absence of treatment
- Air Embolism
- Arrhythmia
- Delay to treatment/therapy
- Device embedded in tissue or plaque
- Exposure to body fluids
- Extravasation
- Foreign body embolism
- Foreign body reaction
- Hematoma
- Hemorrhage/blood loss/bleeding
- Hypersensitivity/allergic reaction
- Hemolysis
- Infections
- Infiltration into tissue
- Lacerations
- Needle stick/puncture
- Nerve damage
- Heparin Induced Thrombocytopenia Pain
- Perforation of vessels
- Phlebitis
- Therapeutic response decreased
- Thrombocytopenia
- Thromboembolism
- Thrombosis/thrombus
- Ulcer
- Unspecified tissue injury
- Vasoconstriction
- Risks Normally Associated with Local or General Anesthesia, Surgery, and Post-Operative Recovery

Kit Contents/Accessories

The primeMIDLINE Catheter kit contains all items to insert the midline catheter into the patient using the modified Seldinger technique.

- MaxPlus Clear Needle-Free Connector 10 ml Syringe
- #11 Safety Scalpel
- Statlock Stabilization Device
- 21 GA x 7 cm Introducer Needle
- Tape Measure
- Tearaway Introducer
- 0.018" Guidewire
- Tuohy Borst Adapter
- Twisted Wire Stylet
- primeMIDLINE Catheter

The primeMIDLINE catheter can be used with fluid delivery accessories via luer connector meeting ISO requirements.

Directions for Use

1. Directions for Modified Seldinger Insertion

1.1. Read instructions carefully before using this device

- Read instructions carefully before using this device.
- The catheter should be inserted, manipulated, and removed by a qualified, licensed physician or other qualified health care professional under the direction of a physician.
- The medical techniques and procedures described in this instructions for use do not represent all medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.
- Use standard hospital protocols when applicable.

MR Information

The primeMIDLINE Catheters are MR-Safe.

2. Prior to Placement

Examine the package carefully before opening to confirm its integrity and that the expiration date has not passed. The catheter is supplied in a sterile package and is non-pyrogenic. Do not use if the package is damaged or opened or if the expiration date has passed. Inspect kit for inclusion of all components. Do not re-sterilize.

☞ **WARNING: Do not use the catheter if there is any evidence of mechanical damage, breakage, or signs of possible leakage. Any damage to the catheter may possibly lead to catheter rupture, fragmentation, air embolism, and surgical removal.**

3. Identify Insertion Site and Vein

Based on patient assessment. Recommended veins are basilic, cephalic and median cubital veins.

☞ **WARNING: (Pediatric) Insertion techniques and placement locations are often modified according to the size and developmental age of the child. Only clinicians experienced in proper positioning and placement of the midlines in pediatric patients should place this catheter in this patient population.**

Ultrasound is recommended for identifying the location of the vein and confirming its diameter. Assess the vasculature without a tourniquet in place.

Assess the selected vein to ensure the vessel size is adequate to accommodate the catheter being placed.

4. Patient Position/Catheter Measurement

- Position the arm at a 90°.
- For midline placement, measure to the desired tip location in the proximal portion of the extremity just distal of the shoulder and deltoid muscle.

5. Skin Preparation

- Don prep gloves.
- Apply underdrape.
- Prepare the site according to institutional policy using sterile technique.

6. Sterile Field Preparation

- Apply tourniquet above the intended insertion site.
- Don sterile gloves.
- Complete sterile field preparation

7. Prepare Catheter

- Strict aseptic technique must be used during insertion, maintenance, and catheter removal procedures. Provide and maintain a sterile field.
- Attach a pre-filled syringe to the luer attachment on the Tuohy Borst adapter
- Pre-flush all catheter lumens with sterile saline to wet the stylet.
- Remove the syringe after pre-flushing.
- Catheters can be trimmed to the desired length if a different length is desired due to considerations of patient size, anatomy and desired point of insertion. Follow hospital protocol.

☞ **WARNING: (Pediatric) This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion, or cardiac tamponade.**

Caution: Catheter tip should terminate at the level of the axilla in children and adults [INS, 2021].

- When trimming the catheter, measure the distance from the zero mark to the desired tip location, catheter markings are in centimeters. Disconnect the stylet assembly and retract the stylet well past the desired point where the catheter is to be trimmed.
- Inspect the cut tip to ensure no fragments or rough edges are present.
- Once the catheter has been trimmed, only advance the stylet to the distal end of the catheter. Never advance the stylet beyond the tip of the catheter. The catheter stylet assembly can now be introduced as described in the following information.

Caution: Never attempt to cut stylet.

Caution: The stylet should be wetted by a flush procedure prior to its withdrawal or repositioning.

Caution: Never close clamp on catheter stylet; stylet and catheter damage may result.

Caution: The stylet should not be advanced beyond the tip of the catheter.

8. Insertion

- Use aseptic technique.
- Anesthetize the insertion site with local anesthesia using facility policies and procedures.
- Insert the introducer needle into the target vein. Observe needle for “flash” to ensure proper placement.
NOTE: See Safecan Safety Introducer Needle instructions for use (below)
- Release tourniquet.
- Insert the flexible end of guidewire into needle and advance guidewire to desired depth into the target vein.

NOTE: See guidewire instructions for use (below)

Caution: The length of the guidewire inserted is determined by the size of the patient and should not be advanced past the axilla.

Caution: Never attempt to cut guidewire.

Caution: Do not insert stiff end of guidewire into vessel as this may result in vessel damage.

- Gently remove **needle**, while holding the guidewire in position.
Caution: If the guidewire must be withdrawn while the needle is inserted, remove both the needle and the guidewire together to prevent the needle from damaging or shearing the guidewire.
- Thread dilator/introducer over the proximal end of the guidewire into target vein. A twisting motion of the dilator/introducer may assist insertion. A small incision may be made adjacent to the guidewire to facilitate insertion of the sheath and dilator. Verify institutional guidelines and policies regarding the use of a safety scalpel.

NOTE: See Tearaway Sheath Introducer instructions for use (below)

NOTE: See safety scalpel instructions for use (below)

Caution: DO NOT bend the dilator/introducer during insertion as bending will cause the sheath to prematurely tear.

- Once the dilator/introducer is in place in the vein, withdraw the dilator and guidewire together as one piece using caution to leave the introducer sheath in position.
Caution: Place a finger over the orifice of the sheath to minimize blood loss and risk of air embolism.
- Insert distal tip of catheter through the sheath until catheter tip is correctly positioned in the target vein.
- Advance the catheter slowly.
- Complete catheter advancement into desired position.

☞ **WARNING: (Pediatric) This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion, or cardiac tamponade.**

Caution: Midline catheters are inserted into a peripheral vein of the upper arm via the basilic, cephalic, or brachial vein with the terminal tip located at the level of the axilla in children and adults. [INS, 2021].

- Stabilize the insertion site by applying gentle pressure to the vein distal to the insertion site.
- Remove the tear-away sheath by slowly pulling it out of the vessel while simultaneously splitting the sheath by grasping the tabs and pulling them apart.

Caution: Do not pull apart the portion of the sheath that remains in the vessel. To avoid vessel damage and skin laceration, pull back the sheath and tear the sheath only a few centimeters at a time. Never leave sheath in place as an indwelling catheter as damage to the vein may occur.

Caution: Do not clamp the catheter using forceps; use only the in-line clamp(s) provided.

- Unscrew the Tuohy Borst adapter from the catheter luer connector.
- Grasp the stylet and Tuohy Borst adapter and remove by slowly pulling back with a constant motion.
Caution: To avoid difficulty and/or bunching of the catheter lumens while removing the stylet, flush catheter prior to removal. The catheter may need to be repositioned to allow for removal of the stylet.
Caution: Do not attempt to reinsert stylet once it has been withdrawn.

Caution: Never leave stylet in place after catheter insertion; injury may occur. Remove stylet after insertion.

- Attach saline filled syringe to lumen(s). Aspirate for adequate blood return and then flush with saline to ensure patency with at least 10 ml of sterile saline to ensure patency. Lock each lumen of the catheter with sterile saline.

Caution: If excessive resistance to blood aspiration is experienced, the catheter may need to be repositioned to obtain adequate flow.

- Attach a new sterile MaxPlus Clear Needle-Free Connector (included in the kit) to each catheter hub.
NOTE: See MaxPlus Clear Needle-Free Connector instructions for use (below)
- Record catheter length and manufacturing lot number on patient's chart.
- Avoid air embolism by keeping catheter tubing clamped at all times when not in use. Aspirate then flush the catheter with saline prior to each use. With each change in tubing connections, purge air from all connecting tubing and accessories.

☞ **WARNING: If signs of extravasation exist, discontinue injections. Begin appropriate medical intervention immediately.**

9. Use of Safecan™ Safety Introducer Needle

Indications for Use:

The Safety Introducer Needle incorporates an active safety clip mechanism to minimize inadvertent needlesticks when needle is used to introduce a guidewire.

Contraindications:

The Safety Introducer Needle should not be used in patients with known hypersensitivity to any of the materials employed.

Risks:

This Safety Introducer Needle is designed to reduce the risk of accidental needlesticks; however, care must be taken to avoid needlesticks. Universal precautions must be adhered to in accordance with Centers for Disease Control and Prevention/Occupational Safety and Health Administration (CDC/OSHA) standards for bloodborne pathogens, when utilizing any needles that penetrates the skin to avoid the risk of exposure to contaminated blood.

Warnings:

- After withdrawal of the Safety introducer Needle from the patient and activation of the safety clip, do not attempt to re-insert the needle into the patient. Also, do not attempt to recap the needle with the protective guard. Dispose of the needle immediately into an appropriate sharps container.
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

Instructions:

NOTE: Both hands should remain behind the needle bevel at all times.

- Inspect prior to use to ensure no damage has occurred during shipping and handling. Do not use if packaging or product is damaged or contaminated or protective caps are loose or missing. Use of damaged or contaminated device might lead to injury, illness, or death of the patient.
- If during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.
- After disinfection of the puncture site and removal of the protective guard from the needle, puncture desired vessel.

- Insert guidewire into the proximal end of the needle hub.
- Advance guidewire through needle into vessel.
- Advancement of guidewire may require gentle rotating motion.
CAUTION: Do not withdraw guidewire against needle bevel to avoid possible severing or damaging of guidewire. If guidewire will not advance, simultaneously remove needle and guidewire from puncture site and restart.
- Advance guidewire in routine fashion to desired depth.
- Hold guidewire in place, withdraw Safety Introducer from the vessel and remove the needle from guidewire.
CAUTION: Do not activate safety clip before removing the needle from the guidewire.
- Grasp clear clip housing, which houses the metal safety clip, with your free hand. Advance clear clip housing forward along the Safety Introducer Needle toward needle bevel.
- In one continuous motion, advance clip housing over needle bevel activating safety mechanism. The clip will remain on the needle bevel. The clip housing will separate from the clip and needle.
- Dispose of needle and clear clip housing immediately into sharps container.

10. Use of .018" x 45cm Galt Medical Guidewire

Indications for Use:

Guidewires associated with this IFU are specifically intended for use in percutaneous procedures, to introduce and position catheters and other interventional devices within the peripheral vasculature.

Contraindications:

Use of the guidewire is contraindicated if the patient has a known or suspected obstruction in the vasculature.

Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

Cautions:

- This procedure should only be performed by physicians thoroughly trained in this procedure.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- If resistance is met when advancing or withdrawing the guidewire, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

Warnings:

- Do not alter this device in any way.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear, unravel, or PTFE coating may scrape off.
- Do not resterilize.

Instructions:

- Peel open package and place contents on sterile field.
- Prep skin and drape in area of anticipated vein-puncture as desired.
- Locate vessel using small gauge needle and syringe (not supplied).
- Insert appropriate gauge Introducer Needle into vessel - observe for flash back.
- Flush the guidewire with sterile heparinized normal saline or a similar isotonic solution. Insert desired flexible tip of the guidewire through the introducer needle into vessel. Slowly advance guidewire to required depth. Advancement of “J” tip may require a gentle rotating motion.
Caution: At no time should guidewire be advanced or withdrawn when resistance is met. Determine cause of resistance before proceeding.
- Hold guidewire in place and remove Introducer Needle.
Caution: Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.
- With the guidewire in place follow the instructions for use provided by the manufacturer of the introducer, catheter or other device, as well as the medical facility's standard procedure for these types of devices.
- To lengthen the flexible segment of a movable core guidewire, gently withdraw the core. To shorten the flexible segment, push the core in carefully.
Caution: When the guidewire is in a vessel, do not advance the movable core if the tip is in a curved shape. Never twist or force the core because excessive force may cause it to penetrate the coil and damage the vessel.

11. Use of Galt Medical Tearaway Sheath Introducer**Indications for Use:**

These introducers are used for the percutaneous introduction of diagnostic or therapeutic devices, such as catheters and pacing leads, into the vasculature.

Contraindications:

Use of the introducer is contraindicated if the patient has a known or suspected obstruction in the vessel. There is increased risk of pneumothorax for the patient who has severe chronic lung disease. Poor healing may result in the patient who has had irradiation to the anterior chest.

Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

Cautions:

- This procedure should only be performed by physicians thoroughly trained in this procedure.
- If resistance is met when advancing or withdrawing the guidewire or the introducer, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- Do not use forceps to break the handle and/or to peel the sheath as this may damage the sheath and cause premature withdrawal of the sheath from the patient.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.

- When locking dilator into the outer sheath, Do Not over-tighten dilator. Damage to the locking feature in the sheath will occur if dilator is rotated more than 1/8 turn in the sheath.
- Symmetrical peeling of the sheath is critical.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

Warnings:

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

Instructions:

- Peel open package and place contents on sterile field.
- Prep skin and drape in area of anticipated vein-puncture as desired.
- Distend the vessel of choice following standard hospital practice for venipuncture. The vein will be much easier to locate if the patient is well hydrated.
- Insert needle into vessel. The needle position should be verified by observing venous blood return.
- The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person.
- Aspirate the puncture needle using the syringe.
- Remove the syringe and insert soft tip of the guidewire through the introducer needle into the vessel. Advance guidewire to required depth. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire position is suggested.
- Hold guidewire in place and remove introducer needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.
- Thread the dilator/sheath assembly over the guidewire.
- Advance the dilator and sheath together with a twisting motion over the guidewire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.
- Remove the vessel dilator and guidewire, leaving the sheath as a conduit into the vessel. Immediately place a finger over the remaining sheath orifice to prevent excessive bleeding or possible air aspiration.
- Advance the catheter or device through the sheath and into the vessel.

Note: If desired, lubricate the device with a light coating of appropriate sterile lubricant to facilitate passage through the introducer sheath. Apply moderate pressure until the device tip pushes through the distal end of the sheath. Some resistance may be felt when advancing a device having tines or a flange tip.

- When the catheter or other device is properly located, peel the sheath away from the device by using the ends of the sheath as handles and pulling them so that the sheath separates longitudinally and breaks at the distal end.

Caution: Do not attempt to withdraw a tined or flanged device through the introducer sheath. If the device must be withdrawn, remove the sheath first and then gently pull the device out of the vessel. Repeat the insertion procedure with a new introducer.

12. Use of Bard Parker Protected Disposable Scalpel

Intended Use:

- Intended for tissue separation and other procedures that require a sharp blade to puncture or cut.

Precautions:

- This device is single use only. Reuse of device could result in infection/contamination and/ or device failure which could lead to patient harm.
- Ruler shown on scalpel handle should be used as an indicator only, and not for accurate measurements.
- During use avoid twisting, bending or putting excessive force or strain on the blade in order to help prevent breakage.
- A singular blade is intended to cut tissue up to 24 inches.
- If blade becomes dull or breaks, dispose and replace product.
- If safety shield or housing malfunction, dispose and replace product.

Instructions:

- To Use
 - Depress top button on protective shield.
 - Retract shield to rear locked position to fully uncover blade.
- To Cover
 - Depress top button on protective shield.
 - Advance shield to rear locked position to fully uncover blade.
- Disposal
 - Dispose in approved puncture-resistance sharps container according to facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

Notice to Users and/or Patients in EU:

- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

13. Use of MaxPlus Clear Needle-Free Connector

Indications for Use:

The MaxPlus is a sterile single patient use positive displacement connector for needleless access to the IV line and/or IV catheter during IV therapy. The MaxPlus connector can be used for direct injection, intermittent infusion, continuous infusion, or aspiration.

Description:

The MaxPlus needleless connector is a closed luer activated device. The accessing ISO male luer from standard administration sets, extension sets, and syringes activate the flow of fluid through the device. The MaxPlus features Tru-Swab® technology which provides a flat, smooth surface for optimum disinfection during pre-access swabbing. The positive displacement feature of the MaxPlus product produces a positive bolus of fluid to clear the catheter upon disconnection of the male luer.

Instructions:

- Without touching the blue top, remove the MaxPlus from the package. Do not contaminate.
- Attach primed IV set or syringe with flush/priming solution to the MaxPlus. Invert the MaxPlus and prime in accordance with facility protocol.
- Remove the protective cover and attach MaxPlus to the hub of the extension set or vascular access device by inserting luer into the hub and rotating until secure. Do not over tighten.
- Prior to every access, always swab the top of the MaxPlus with an appropriate antiseptic and allow to dry.
- Attach luer from primed IV set or syringe to the MaxPlus. If luer is a two piece spin collar, first pull back the collar and insert luer in a straight in motion and rotate ¼ turn clockwise; then push the spin collar forward and tighten. If syringe has a luer slip, insert and rotate ¼ turn clockwise to secure connection. Do not leave luer slip unattended.
- Flush the MaxPlus after each use with normal saline or in accordance with facility protocol.
- To disconnect male luer from MaxPlus, hold MaxPlus and rotate mating luer until disconnected. Upon removal of the luer from the MaxPlus a positive displacement of fluid will occur to reduce retrograde flow into the catheter lumen. If blood or fluid residuals remain on the connector surface after access, swab surface.
- Apply clamp, after disconnecting luer.

Warnings and Cautions:

- The device should not be used with needles, blunt cannula systems, non ISO luer connections, or luer connections with visible defects. Doing so may result in leakage and/or failure of the device.
- Luer slip connections should not be left unattended due to potential for disconnection.
- For proper use, clinicians must be familiar with and trained in the use of the MaxPlus device. Its use should be preceded by an established facility protocol. The device should be disinfected with an appropriate antiseptic agent, such as 70% IPA, prior to each access.
- The device is intended for use with ISO luer lock and luer slip connectors provided on standard IV administration sets, extension sets, and syringes.
- Trace lines before connection. Verify the line being connected to the MaxPlus is the appropriate Intravenous Therapy line.
- Failure to properly prime the device can result in reflux.
- The MaxPlus should be changed according to facility protocol or in accordance with currently recognized guidelines for IV therapy.

- Sterile disposable single patient use device may be accessed multiple times per hospital protocol. Reuse, reprocessing or re-sterilizing may lead to patient infection or other illness/injury.
- Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.

14. Catheter Securement and Wound Dressing:

- Secure the catheter using the provided StatLock® stabilization device and/or per institutional policy.
NOTE: See StatLock® instructions for use (below)
- Dress the catheter with an occlusive dressing per institutional policy.

15. Use of StatLock®:

Device Description

- The StatLock™ Stabilization Device can be used with various marketed catheters. The device secures the catheter or medical tube to minimize the risk of migration and dislodgement. The device may be used up to 7 days. Natural rubber latex is not part of the material formulation.

StatLock® Intended Use

- The StatLock™ stabilization device is a stabilization device for compatible medical tubes or catheters.

StatLock® Intended User

- Procedure must be performed by an experienced healthcare professional with knowledge of anatomical landmarks, safe technique and potential complications.

StatLock® Intended Patient Population

- The StatLock™ stabilization device is intended for patients who need catheter or other medical tubing securement.

StatLock® Clinical Benefit

- The StatLock™ Stabilization Device provides stabilization of the catheter to the skin to minimize migration and dislodgement of inserted catheters.

StatLock® Contraindications

- Known tape or adhesive allergies

StatLock® Warnings

- Do not use the StatLock® stabilization device where loss of adherence could occur, such as with a confused patient, diaphoretic or nonadherent skin, or when the access device is not monitored daily.
- This is a single use device. Re-use and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.
- After initial sterilization do not resterilize. The sterility of the single use device is not guaranteed following resterilization because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Resterilization may compromise the structural or design characteristics and may lead to an unpredictable loss of functionality and/or device failure.
- May cause anaphylactic shock if patient is allergic to one or more of the materials.
- Use may induce skin maceration or damage; improper removal may cause patient discomfort.
- Use of device could lead to compression, kinking, or pinching of the catheter/tubing. May cause reduced flow through the catheter and may lead to hemolysis.
- Non-sterile product use may lead to infection.

- Do not alter the StatLock™ stabilization device or components.

StatLock® Precautions

- Please contact BD to obtain sterilization process information.
- Observe universal blood and body fluid precautions and infection control procedures, during application and removal of the StatLock® stabilization device.
- Remove oil and moisturizer from targeted skin area.
- Avoid contact with alcohol or acetone. Both can weaken bonding of components and the stabilization device pad adherence.
- Minimize catheter/tube manipulation during application and removal.
- Monitor daily and replace when clinically indicated, at least every 7 days. Use aseptic technique for application and removal.
- Care should be taken in removing catheter from retainer to prevent inadvertent access line dislodgement or tip malposition; may lead to hematoma or bleeding.
- Leakage of fluid may cause skin irritation.

Directions for use

Use aseptic technique for application and removal of the StatLock™ Stabilization Device

Application:

- Prepare
 - When included in the kit, place adhesive foam strip over catheter hub but NOT over the insertion site. A secondary foam strip may be used to fashion a J-loop with the extension set/IV tubing once the rest of the StatLock™ application is completed.
 - Cleanse and degrease insertion and stabilization site with facility approved anti-microbial agent and allow to dry completely.
 - Then apply included skin protectant preparation pad solution for enhanced adherence and skin protection. **Allow to dry completely.**
- Press
 - Where applicable, align anchor pad so directional arrow(s) points toward the insertion site.
 - Press the catheter or extension set into the StatLock™ stabilization device retainer, while providing support underneath.
 - Ensure the retainer is firmly attached to the fitting and will not allow the catheter or extension set to slip.
 - Where applicable, close the retainer and verify that the catheter is secure in the stabilization device.
- Peel
 - Peel off the StatLock™ stabilization device paper backing or liner.
- Place
 - Place the StatLock™ stabilization device on targeted skin that has been properly prepared.

Removal:

Where applicable, carefully remove overlying dressing using manufacturer's instruction or facility policy.

- Dissolve
 - Dissolve adhesive with alcohol swab while gently lifting the StatLock™ stabilization device pad.
- Disengage
 - Disengage catheter/tube from the StatLock™ stabilization device retainer carefully, to avoid inadvertent line dislodgment.

- Document
 - Document the StatLock™ stabilization device dressing change in the patient chart.
- Dispose.
 - Dispose of all equipment in appropriate containers.

Note: If placing a new StatLock™ stabilization device, the foam strip should be replaced after the main device change is completed. If the catheter is being taken out, the foam strip should be removed just before catheter removal.

Note: Users should report any serious incident related to the device to the Manufacturer and National Competent Authority.

16. Power Injection Procedure

- Remove the injection/needleless cap from the *prime*MIDLINE Catheter, if present.
- Using a 10 ml or larger syringe(s), aspirate for blood return and flush catheter lumen to assure free flowing patency. Discard syringe(s).
- Detach syringe

☞ **WARNING: Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure.**

- Contrast media should be warmed to body temperature prior to power injection.

Power Injection Flowrate Test Details	
Extension tube length and internal diameter	Length: 2.738 in ± 0.125 in
	ID: 0.066 in ± 0.002 in
Injectate dynamic viscosity	11.8 mPa·s
Injectate temperature	37°C

- Attach the power injection device to the *prime*MIDLINE Catheter per power injection device manufacturer's recommendations.
- Complete power injection study taking care not to exceed the flow rate limits below:

Catheter Size	Maximum Flow Rate	Maximum Pressure
3F (1.00 mm) Single Lumen	0.36 L/min (6 ml/s)	2068 KPa (300 psig)
4F (1.33 mm) Single Lumen	0.42 L/min (7 ml/s)	2068 KPa (300 psig)
5F (1.67 mm) Single Lumen	0.42 L/min (7 ml/s)	2068 KPa (300 psig)
5F (1.67 mm) Dual Lumen	0.42 L/min (7 ml/s)	2068 KPa (300 psig)

☞ **WARNING: The pressure limiting feature of the power injector may not prevent over-pressurization of an occluded catheter.**

☞ **WARNING: Exceeding the maximum indicated flow rate or maximum pressure of 300psi may result in catheter failure.**

- Disconnect the power injection device.
- Flush the *prime*MIDLINE Catheter with 10 ml of sterile normal saline, using a 10 ml or larger syringe. For multi-lumen catheters, flush all lumens after power injection.
- Replace the injection/needleless cap on the *prime*MIDLINE Catheter.

17. Infusion

- All connections should be examined carefully before infusion begins.
- Disinfect needleless connector.
- Unclamp extension leg clamps.
- Confirm catheter patency.
- Connect infusion set/syringe.
- Complete infusion as prescribed.
- Disconnect infusion set/syringe.
- Disinfect needleless connector.
- Connect 10ml syringe or larger filled with sterile saline.
- Confirm catheter patency.
- Clamp extension leg clamps.
- Disinfect needleless connector.

18. Blood Sampling

- Disinfect needleless connector per institutional policy if connector is indicated for blood sampling.
- Disconnect needleless connector and disinfect catheter luer hub if connector is not indicated for blood sampling.
- Confirm catheter patency per institutional policy.
- Attach blood withdrawal device to either needleless connector or catheter luer hub.
- Unclamp extension leg clamps on catheter.
- Perform blood sampling per institutional policy.
- Clamp extension leg clamps on the catheter.
- Disconnect blood withdrawal device.
- Confirm catheter patency per institutional policy.
- Disinfect catheter luer hub or needleless connector.
- Attach a new needle free connector per institutional policy, as needed.

19. Catheter Maintenance

- Dressing Changes - Assess the dressing in the first 24 hours for accumulation of blood, fluid or moisture beneath the dressing. During all dressing changes, assess the external length of the catheter to determine if migration of the catheter has occurred. Periodically confirm catheter placement, tip location, patency and security of dressings. An occlusive dressing should cover the exit site at all times. The dressing should be changed per institutional policy or any time the dressing becomes soiled, wet, or non-occlusive.

NOTE: When using alcohol or alcohol containing antiseptics with the *prime*MIDLINE Catheter, care should be taken to minimize contact duration. Solutions should be allowed to completely dry before applying an occlusive dressing. Chlorhexidine gluconate is the suggested antiseptic to use. Follow hospital protocol.

☞ **WARNING: Minimize exposure of catheter to alcohol.**

☞ **WARNING: Acetone and polyethylene glycol containing ointments should not be used with the *prime*MIDLINE Catheter, as these may cause failure of the device.**

☞ **WARNING: Do not use scissors to remove dressings to minimize the risk of cutting the catheter.**

20. Flushing

- Follow institutional policy for flushing frequency.
- The catheter should be flushed with normal saline prior to drug administration.
- After drug administration, each lumen should be flushed again with normal saline.
- Flush with 20 ml sterile saline after blood therapy.
- Injection cap(s) or needleless access port(s) should be changed per institutional policy.

☞ **WARNING: If signs of extravasation exist, discontinue injections. Begin appropriate medical intervention immediately.**

21. Catheter Performance

- If resistance is encountered when aspirating or flushing, the lumen may be partially or completely occluded.

☞ **WARNING: Do not flush against resistance.**

- If the lumen will not aspirate or flush, follow institutional guidelines to determine if catheter is occluded with blood or precipitate and follow appropriate recommended procedure.

Cleaning Catheter Exit Site

- Follow institutional guidelines and policy.

22. Catheter Removal

☞ **WARNING: Only a clinician familiar with the appropriate techniques should attempt the following procedures.**

Caution: Always review facility protocol, potential complications and facility treatment, warnings, and precautions prior to catheter removal.

- Remove dressing and catheter stabilization device.

Caution: Do not use scissors to remove dressing to minimize the risk of cutting the catheter.

- Grasp catheter near insertion site.
- Remove slowly. Do not use excessive force.
- If resistance is felt, stop removal. Apply a warm compress and wait 20-30 minutes before resuming removal procedure.
- Inspect the end of the catheter and the catheter length to ensure complete removal.

Caution: After use, this product may be a potential biohazard. Handle and dispose of catheter in accordance with medical practice and applicable local, state, and federal laws and regulations.

☞ **WARNING: Catheter is intended for single use. DO NOT REUSE. Reuse and/or repackaging may create a risk of patient / user infection or compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure and/or injury, illness, or death of the patient.**

23. Disposal After Use

- After use, medical products and accessories pose a potential biological hazard. For this reason, the products and their accessories should be handled and disposed of in accordance with recognized medical procedure, and in compliance with applicable local, state, and federal laws and regulations.

Notice

Report any serious incident that has occurred in relation to the device to PFM Medical, Inc. and the competent authority of the Member State in which the user and/or patient is established.

Warranty

PFM Medical, Inc. warrants that this product was manufactured according to applicable standards and specifications. Patient condition, clinical treatment, and product maintenance may affect the performance of this product. Use of this product should be in accordance with the instructions provided and as directed by the prescribing physician.

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