

pfmmedical
Quality and Experience

ASEPT[®] Select

ASEPT[®] Select Drainage Kit

Offering patients the benefit
of both gravity and vacuum
drainage.


www.pfmmedicalusa.com



The New Standard in Out-Patient Drainage

The ASEPT Select Drainage Kit offers patients the benefit of both gravity and vacuum drainage. The vacuum provides the patient efficient and effective fluid drainage, while gravity helps the patient have a more pleasant drainage experience. Color coded slide clamps allow the patient to switch from vacuum to gravity during the procedure at any time. The ASEPT Select Drainage Kit can be used for both chest and abdomen drainage in the out-patient setting.

Features



- 1 Shatter-proof plastic bottle
- 2 Two color coded slide clamps to choose between gravity and vacuum
- 3 Disk check valve to help prevent potential air and fluid back-flow
- 4 Luer-lock connection to ASEPT catheter
- 5 Easy to use multi-step pinch clamp to conveniently regulate fluid flow

Procedure Pack

Each ASEPT Select Drainage Kit comes with an ASEPT Procedure Pack, which contains all the items necessary to perform an aseptic drainage procedure. The ASEPT Procedure Pack includes:

- Foam Catheter pad
- 4x4 Gauze Pads
- 3 Alcohol swabs
- 2 Pairs of Gloves
- Self-Adhesive Surgical Dressing
- Emergency Clamp



Ordering Information

Item No.	Description	Case Qty
M7006	ASEPT Select 600mL Drainage Kit (Includes low vacuum bottle and Procedure Pack)	10/Case
M7007	ASEPT Select 1000mL Drainage Kit (Includes low vacuum bottle and Procedure Pack)	10/Case

www.pfmmedicalusa.com

pfmmedical

pfm medical, inc.

1916 Palomar Oaks Way, Suite 150 • Carlsbad, CA 92008

Ph: 760.758.8749 • Fx: 760.758.1167

customerservice@pfmmedicalusa.com • www.pfmmedicalusa.com

ASEPT is a trademark of pfm medical, inc. Registered in U.S. Patent and Trademark Office