

**primePICC**  
**Product and Kit Configuration**

Component	Basic Kit	Nursing Kit (N Suffix)	Long Wire Kit (L Suffix)
Catheter length: Lumen Configurations: S = Single   D = Dual   T = Triple • all catheters are CT rated at 5 ml/sec and 300 psi	3F: 45 cm   Lumen: S	3F: 45 cm   Lumen: S	3F: 45 cm   Lumen: S
	4F: 50 cm   Lumen: S	4F: 50 cm   Lumen: S	4F: 50 cm   Lumen: S
	5F: 55 cm   Lumen: S, D, T	5F: 55 cm   Lumen: S, D, T	5F: 55 cm   Lumen: S, D, T
	6F: 60 cm   Lumen: D, T	6F: 60 cm   Lumen: D, T	6F: 60 cm   Lumen: D, T
Tuohy Borst Adapter w pre-loaded .015" x 75 cm stylet	Yes	Yes	No stylet
21G x 7 cm Introducer Needle w Echo Tip	Yes	No	Yes
21G x 7 cm Safety Needle w Echo Tip	No	Yes	No
Guidewire	.018" x 80 cm Nitinol w Gold Tip	.018" x 45 cm Nitinol w Tungsten Tip	.018" x 130 cm Nitinol w Gold Tip
Tearaway Introducer	3F: 3.5F x 5 cm	3F: 3.5F x 5 cm	3F: 3.5F x 10 cm
	4F: 4.5F x 10 cm	4F: 4.5F x 5 cm	4F: 4.5F x 10 cm
	5F: 5.5F x 10 cm	5F: 5.5F x 5 cm	5F: 5.5F x 10 cm
	6F: 6.5F x 10 cm	6F: 6.5F x 5 cm	6F: 6.5F x 10 cm
Neutral Displacement Needlefree Connector (1 per lumen)	Yes	Yes	Yes
Safety Scalpel	Yes	Yes	Yes
Statlock® PICC Plus Securement Device	Yes	Yes	Yes
Tape Measure	Yes	Yes	Yes