


Veta™ Peritoneal Dialysis Catheter Kit

Accurate placement
through visualization
with a minimal incision

www.pfmmedical.com



Features and Benefits



- ▶ **Silicone Material with Radiopaque Stripe**
 - ▶ Minimizes skin irritation at exit site
 - ▶ Used for long term implantation
 - ▶ Helps pinpoint catheter location
 - ▶ Large spiral flow holes




- ▶ **Catheter Cuffs**
 - ▶ Provide tissue in-growth for anchoring the catheters
 - ▶ Provide a barrier against possible infection
 - ▶ Prevent leakage of dialysate along the catheter track



- ▶ **Introducer Kit**
 - ▶ Safety scalpel
 - ▶ Guide wire
 - ▶ Introducer and needle
 - ▶ Catheter introducer and dialator
 - ▶ Syringe
 - ▶ Tunneler
 - ▶ Gauze pads
 - ▶ Red flow clamp
 - ▶ Catheter hub and end cap
 - ▶ Catheter curled 62cm

Accessories



- ▶ **Stylet**
 - ▶ Malleable stainless steel stylet
 - ▶ Straightens catheter for ease of placement

Intended Use

- ▶ **Acute and Chronic Access to the Peritoneal Cavity**

Ordering Information

Veta™ Peritoneal Dialysis Catheter

Item No.	Description	Length	Box Qty
M7100	Straight, Double Cuff	42cm	5/Box
M7200	Curled, Double Cuff	62cm	5/Box
M7225	Kit with Catheter Curled	62cm	5/Box
M7230	Kit (Only)	-	5/Box
M8001NS	Catheter Stylet (Non-Sterile)	-	1/Box
M8002	Catheter Stylet (Sterile)	-	5/Box

www.pfmmedical.com

pfmmedical

pfm medical, inc.
 1815 Aston Ave., Ste. 106 • Carlsbad, CA 92008
 Ph: 760.758.8749 • Fx: 760.758.1167
customerservice@pfmmedicalusa.com • www.pfmmedical.com