

pfmmedical

Quality and Experience

Veta™ Peritoneal Dialysis Catheter

**For Acute and Chronic
Access to the Peritoneal
Cavity**


www.pfmmedical.com



Features and Benefits



- ▶ **Silicone Material**
 - ▶ minimizes skin irritation at exit site
 - ▶ used for long term implantation



- ▶ **Radiopaque Stripe**
 - ▶ helps pinpoint catheter location



- ▶ **Large Spiral Flow Holes**
 - ▶ maximize fluid delivery and drainage




- ▶ **Coiled Distal Catheter End**
 - ▶ helps catheter maintain position in pelvic area

- ▶ **Catheter Cuffs**
 - ▶ provide tissue in-growth for anchoring the catheters
 - ▶ provide a barrier against possible infection
 - ▶ prevent leakage of dialysate along the catheter track

Intended Use

- ▶ **Acute and Chronic Access to the Peritoneal Cavity**

Ordering Information



Veta™ Peritoneal Dialysis Catheter

Item No.	Description	Length	Box Qty
M7100	Straight, Double Cuff	42 cm	5/Box
M7200	Curled, Double Cuff	62 cm	5/Box

www.pfmmedical.com